



## Simply Easy to Use.

With 30 years of experience in lifting and transferring people in care, Liko™ has applied its proven technology and design excellence to create the new Liko™ M220 and M230 mobile lifts.

These easy-to-use, innovative lifts are designed to meet standard lifting needs in home care and nursing home environments with the same durability and quality you've come to trust from Liko™.

Like many other Liko™ products, the Liko™ M220 and M230 mobile lifts feature an intuitive design that makes the lifts extraordinarily easy to use, even for caregivers new to lift operation. Coupled with a wide range of available slings and accessories, the lifts create a secure and comfortable lifting experience for both residents and caregivers.

The Liko™ M230 model is fully electric, including the base-width adjustment, while the Liko™ M220 model has a manual base.

### Important Features

- Easy-to-use, intuitive controls
- Long lift arm provides full access to center of bed
- Easy to clean
- Electrical lifting and base-width adjustment (manual base-width adjustment on M220 model)
- 182 kg (400 lbs) lifting capacity
- Compatible with world's most complete range of slings enabling you to adapt the lift to the individual resident's particular needs
- Can be used with a scale for weighing (accessory)
- Attractive design and competitive price

# Liko™ M220 and M230 Mobile Lifts

## The Liko™ M220 and M230 mobile lifts at a glance



*With the Liko™ M220 and M230 mobile lifts, you can meet all of the most common lifting needs, including lifting from the floor.*



*The long lift arm allows for easy patient rotation and easily lifting the patient to the center of the bed.*

- Lifting capacity: 182 kg (400 lbs)
- Lifting interval: 570-1720 mm (22.4-67.7 inch.)
- Wheels: Front: 75 mm (3 inch) twin wheels  
Rear: 75 mm (3 in) twin wheels with brakes
- Material: Steel
- Emergency lowering: Manual and Electrical
- Models: Liko™ M220 (manual base)  
Liko™ M230 (electric base)

Liko™ M220 and M230 have been tested by an accredited testing institute and comply with the requirements for Medical Device Directive Class I products (MDD 93/42/EEC). Liko™ M220 and M230 comply with the requirements in EN ISO 10535, IEC 60601-1, IEC 60601-1-2, IEC 60601-1-11, ANSI/AAMI ES60601-1, and CAN/CSA C22.2 no 60601-1. Design and Quality by Liko in Sweden. Liko™ is quality certified in accordance with ISO 9001 and its equivalent for the medical device industry, ISO 13485. Liko™ is certified in accordance with environmental standard ISO 14001.



Hill-Rom and Design and Hill-Rom word mark are trademarks or registered trademarks of Hill-Rom Services Inc. Liko and Design and Liko word mark are trademarks or registered trademarks of Liko R&D AB.



Enhancing Outcomes for Patients and Their Caregivers.™

France \_\_\_\_\_ +33 (0)2 97 50 92 12  
 United Kingdom \_\_\_\_\_ +44 (0)1530 411000  
 Deutschland \_\_\_\_\_ 0800 8 63 55 37  
 Nederland \_\_\_\_\_ +31 (0)347 32 35 32  
 Italia \_\_\_\_\_ + 39 (0)2 950541  
 Suisse/Schweiz \_\_\_\_\_ +41 (0)21 706 21 30  
 (deutschsprachig) \_\_\_\_\_ 0800 8 63 55 37

Österreich \_\_\_\_\_ 0800 8 63 55 37  
 Ireland \_\_\_\_\_ +353 (0)1 413 6005  
 Iberia \_\_\_\_\_ +34 (0)93 6856000  
 Nordic Countries \_\_\_\_\_ +46 (0)20 78 10 30  
 Export \_\_\_\_\_ Contact your local Hill-Rom distributor or contact your Area Manager via website or call \_\_\_\_\_ +1 812 934 8173

[www.hill-rom.com](http://www.hill-rom.com)  
[www.liko.com](http://www.liko.com)

Hill-Rom® reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom® makes is the express written warranty extended on the sale or rental of its products.  
 © 2013 Hill-Rom® SARL, ALL RIGHTS RESERVED.