

Procedural Stretcher - P8000

Smoother transport for patients and caregivers



Features and benefits

Steering Plus™ directional system incorporates fully-locking five castor braking, providing improved steering, stopping, cornering and control

An **integrated IV pole** cradle and ergonomic push handles facilitate easier patient transport

Fully integrated siderails mean safer patient transport

Power cord hook for easier storage and reduced caregiver strain during patient transport

Enhancing outcomes for
patients and their caregivers:

Hill-Rom

Procedural stretcher - P8000



Power cord hook



Optional: 4 brake & steer pedals



Optional: push handles and integrated IV pole



Optional: integrated handbrake system

Technical Specifications (standard configuration)

Safe Working Load	317.5 kg
Weight (without surface/accessories)	120.2 kg
Height adjustment	52.1 to 87 cm
Trendelenburg/Reverse Trendelenburg	- 18 to + 18 °
Protection class	IPX4
Deck width	
standard	66 cm
optional	76.2 cm
Width	
siderails stored	74.3 or 84.5 cm
siderails up	81.3 or 91.4 cm
Deck length	190.5 cm
Length	210.8 cm
Surface height	
standard	7.6 cm
optional	10.2 or 12.7 cm

Siderails

height	36.8 cm
length	119.4 cm

Clearance

2.9 to 8.9 cm

Castors

4 x 20.3 cm
single band

Additional options

- Push handles
- Integrated handbrake system
- Hydraulic control pedals on both sides at foot end
- 4 brake & steer pedals
- Integrated O² tank holder / utility tray
- Radiolucent surface
- Upright cassette holder
- Integrated IV pole - 3 stage
- Knee gatch
- Backsaver Fowler™ (includes knee gatch)
- Auto Contour™ (includes knee gatch and Backsaver Fowler)

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Medical devices shown in this brochure are intended for use with patients in departments of healthcare establishments.

These medical devices are regulated health products which, pursuant to such regulation, bear a CE mark. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN145302-02, 18 March 2016

Enhancing outcomes for patients and their caregivers:

Hill-Rom