

Overbed table – TA270

Sturdy functionality with ergonomic design



TA270 TT (tilting function)



TA270 FT (fixed table)

Features and benefits

- The Hill-Rom Autolift™ system offers optimal patient/resident and caregiver safety by reducing the risk of entrapment and overturning – with a simple push, the table top will lift and move securely into place
- Efficient H-shaped undercarriage means the table is compatible with all care environments, delivering maximum working stability
- Ergonomic design features a reduced-width front wheelbase, allowing it to be positioned in all areas of the patient/resident environment
- Locking dual-directional tilted position offers patient/resident an ideal reading position
- The fully-rotational table top and epoxy-coated undercarriage facilitate easy cleaning

Overbed table – TA270

Technical Specifications **model TA270**

Height adjustment (without rails) 730 – 1,080 mm

Width 960 mm

Depth 490 mm

Surface area

model TA270 TT
tilting portion:
740 x 450 mm
fixed portion:
220 x 450 mm

model TA270 FT
fixed table:
960 x 450 mm

Weight 18.5 kg

Maximum load 16 kg

Material

surface area High Pressure
Laminate (HPL)
undercarriage Metallic grey epoxy
coated steel with
plastic cover

Castors

4 x 50 mm
double-band
with 2 brakes

Tilting range

model TA270 TT
max. 30° locked
positions /360° full
rotational

Colours

see colours sheet

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries.
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN150324-01 (revision 2), 07 June 2017

Enhancing outcomes for
patients and their caregivers:

Hill-Rom