

Bedside cabinet with fridge- CH242C1

Offering innovative design for today's patient environment



model shown with optional 4 castors

Features and benefits

- Suitable across all budgets with a wide variety of colour options
- Long-lasting, durable construction, using materials that are both easy to clean and maintain
- 30 litre built-in refrigerator is ideal for keeping chilled refreshments close at hand
- Refrigerator door (along with cabinet door) can be configured to open bilaterally

Bedside cabinet with fridge– CH242C1

Technical Specifications **model CH242C1**

Height	875 mm
Width	510 mm
Depth	460 mm
Weight	~45 kg
Material	
top plate	laminate
body	melamine
Castors	
standard	2 x 50 mm with 2 end fittings on simple base
optional	4 feet in solid beech with 2 x 50 mm on back
optional	4 x 75 mm double-band with 2 brakes
Configuration	fridge (30 litre) 1 drawer and 1 door
Colours and handles	see colours and handles sheet
Additional options	Lockable drawer

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health**. Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers**.

Not all products/options are available in all countries.
For further information about our products and services, please
contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN150369-01, 20 April 2016

Enhancing outcomes for
patients and their caregivers:

Hill-Rom