

Overbed table – TA529

Sturdy functionality with ergonomic design



TA529 (tilting function)

Features and benefits

- Efficient U-shaped undercarriage means the table is compatible with all care environments, delivering maximum working stability
- Ergonomic design features a reduced-width front wheelbase, allowing it to be positioned in all areas of the patient/resident environment
- Epoxy-coated undercarriage facilitates easy cleaning

Overbed table – TA529

Technical Specifications		model TA529 F	model TA529 G	model TA529 K	model TA529 L
Height adjustment (without rails)		700 – 1,000 mm	700 – 1,000 mm	700 – 1,000 mm	700 – 1,000 mm
Width		890 mm	894 mm	890 mm	890 mm
Depth		468 mm	460mm	450 mm	450 mm
Surface area	fixed table:	890 x 468 mm	-	890 x 450 mm	-
	tilting portion:	-	675 x 460 mm	-	670 x 450 mm
	fixed portion:	-	219 x 460 mm	-	220 x 450 mm
Weight		15 kg	15 kg	15 kg	15 kg
Maximum load		15 kg	15 kg	15 kg	15 kg
Material	surface area:	laminate	laminate	High Pressure Laminate (HPL)	High Pressure Laminate (HPL)
	undercarriage:	welded epoxy-coated steel			
Castors		4 x 50 mm double-band with 2 brakes			
Tilting range		-	max. 30° locked positions	-	max. 30° locked positions
Colour undercarriage		light grey			
Colours		see colours sheet			

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health**. Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers**.

Not all products/options are available in all countries.
For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:
Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN150331-02 (revision 2), 06 June 2017

Enhancing outcomes for
patients and their caregivers:

Hill-Rom