

Viking™ L mobile lift

For additional strength



With a lifting capacity of 250 kg, Viking L is ideal for most common lifting needs, and offers a flexible lifting solution. The extensive application area means that for lifting medium-weight or larger patients, the Viking L is suitable in most situations, e.g. lifting from the floor, to/from the toilet, horizontal lifting, and/or gait training.

Features and benefits

There are **dual control options** for operating the lift - either hand control or control panel

Optional quick release hook system means caregivers can easily and rapidly change a sling bar, (add a Liko Stretch for horizontal lifting) or install a scale

Features a **lightweight, environmentally-friendly NiMH battery**

Optional armrests help **support patients during ambulation** as well as making it easier for caregivers to manoeuvre the lift

An extensive assortment of accessories is available including slings, slingbars and scales

For **additional flexibility** there is an optional wireless hand control available

Viking L mobile lift



Handles



Operational panel



Optional Universal TwinBar 670 slingbar



Wired or wireless hand control

Technical Specifications

| | |
|-------------------------------------|-------------------------------|
| Maximum Load | 250 kg |
| Weight | |
| with wired handcontrol | 37.4 kg |
| with wireless handcontrol | 36.5 kg |
| Height | 1390-2030 mm |
| Protection class | IPX4 |
| Slingbar options | |
| standard | Universal Slingbar 450 |
| optional | Universal TwinBar 670 |
| Base | |
| height | 165 mm |
| length | 1360 mm |
| width | 695 - 1070 mm |
| expansion | electric |
| Castors | |
| front | 100 mm double band |
| rear | 125 mm double band |
| Material | aluminium |
| Emergency lifting | electric |
| Emergency lowering | manual and electric |
| Lifting speed (without load) | 21 mm / sec or 13 mm / sec |
| Lifting interval | 1330 mm |

| | |
|-------------------------------|-------------------|
| Maximum lifting height | 1790 mm |
| Sound level | 51 dBA |
| Battery | NiMH |
| Charger | built in |
| Hand control | wired or wireless |
| Product number | |
| with wired hand control | 2040034 |
| with wireless hand control | 2040014 |

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Viking L

Viking is a versatile lift intended mainly for use in health care, intensive care and rehabilitation. Class I
Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2016 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.
Doc. No: 5EN183308-01, 20 July 2016

Enhancing outcomes for patients and their caregivers:

Hill-Rom