

Bedside cabinet – high

Contemporary bedside cabinets that offer versatility and durability for today's patient environment.



B01



B02



B03



B04

Features and benefits

- A varied selection of drawer and door configurations, suitable for all budgets
- Soft-close, shock-absorbent drawers that glide effortlessly into place
- Long-lasting, durable construction, using materials that are both easy to clean and maintain
- Aptly sized to suit the needs of today's patient environments

Bedside cabinet – high


Catoni Collection - Classic



Optional: Sliding tray



Optional: Accessory bar holder



Optional: Simple base



Optional: Two half doors



Optional: Fixed strip with optional colour contrast



Optional: Colour contrast for fixed shelf

Technical Specifications model CH700B4

Height	885 mm
Width	401 mm
Depth	401 mm
Weight	~22 – 29,5 kg
Material	
Body	Melamine
Top plate	High Pressure Laminate (HPL)
Hinge for B01, B02, B03	
standard	on the left side
optional	on the right side
Fixed strip	
standard	for B02, B03, B04
optional	for B01

Castors

standard	4 x 50 mm with brakes
optional	simple base (2 x 50 mm with brakes and 2 end fittings)

Handles/Colours

see handles and colours sheet

Additional options

- sliding tray
- accessory bar holder for one or two sides
- two half doors for B01, B02, B03
- interior removable shelf for B01, B03
- contrasting colours for
 - sliding tray and/or fixed shelf
 - fixed strip and/or fixed shelf

Hill-Rom is a leading global medical technology company with more than 7,000 employees in over 100 countries. We partner with health care providers by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Clinical Workflow, Surgical Safety and Efficiency, and Respiratory Health.** Hill-Rom people, programs, and product brands work towards one mission: **Every day, around the world, we enhance outcomes for our patients and their caregivers.**

Not all products/options are available in all countries. For further information about this product or a service, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

This document is destined solely for use by healthcare professionals. Devices shown in this brochure are intended for use with patients in departments of healthcare establishments. Hill-Rom recommends you carefully read the detailed instructions for safe and proper use included in the documents accompanying the devices. The personnel of the healthcare establishments are responsible for the proper use and maintenance of these devices.

Furniture manufacturer:

Hill-Rom SAS, B.P. 14 - Z.I du Talhouët, 56330 Pluvigner - France (also European Representative)

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2015 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN150313-01, 18 June 2015

Enhancing outcomes for patients and their caregivers:

Hill-Rom